

Bexen cardio

Ready for life

Reanibex 300 AED

Automated External Defibrillator



2 versions:

1. Icons are displayed
2. Icons or ECG are displayed

Biphasic up to 360J

Universal AED algorithm for adult and paediatric patients

CPR guidance device in real time (Push Pad)

Possibility to configure the screen on high contrast mode

“Autotraining” mode and “trainer” option with remote control

Reanibex 300 AED.

User friendly and intuitive.



FEATURES

- Biphasic technology.
- Universal AED algorithm for adult and paediatric patients.
- Large colour display 5.7".
- Possibility of setting the screen on high contrast mode.
- "Push Pad" device, reports the quality of the frequency and depth of chest compressions in real time (optional).
- Autotraining mode (optional).
- Trainer option with remote control (optional).
- Configurable CPR protocol.
- Configured according to ERC/ AHA guidelines 2015 and its review of 2017.
- Automatic daily selftest ensures the device continual readiness.
- USB connection.
- Rugged and durable design for any environment.

ACCESSORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Carrying case (1).
- Ambulance bracket (EN 1789:2007 certified) (2).
- "Push Pad" device.
- Rechargeable and non rechargeable battery (3).
- External battery charger.
- External dual battery charger 12V (DC) or 220V(AC).
- "REANIBEX Data Manager" Software application to manage and organize all the data gathered.
- Rechargeable battery for autotraining and trainer mode (4).



Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobuelta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com

