

Bexen cardio  
Ready for life



## Reanibex 500 EMS

The most compact  
multiparameter  
Monitor  
Defibrillator of  
the market

Biphasic Manual and AED  
up to 360J

Universal AED algorithm for  
adult and paediatric patients

Non invasive pacing

SpO2

EtCO2

12 leads

Up to six ECG waveforms

CPR guidance device in real  
time ( Push Pad)

Data transmission  
(GSM, Bluetooth, USB...)

# Reanibex 500 EMS. Compact. Lightweight.



## FEATURES

- Biphasic technology.
- Universal AED algorithm for adult and paediatric patients.
- Lightweight unit with multiparameter monitoring capabilities.
- 4+6 leads patient cable.
- Synchronized Cardioversion.
- Non-invasive pacing (optional).
- Pulse Oxymetry with Masimo SET for adult and paediatric patients (optional).
- Capnography from Masimo SET (optional).
- Respiration rate (optional).
- “Push Pad” device, reports the quality of the frequency and depth of chest compressions in real time (optional).
- Configured according to ERC/ AHA guidelines 2015 and its review of 2017.
- Automatic daily selftest ensures the device continual readiness.
- Data transmission (GSM, Bluetooth, USB) to EMS, Hospitals... configured according to the user needs (optional).
- Rugged and durable design for any environment.

## ACCESSORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Adult and paediatric ECG electrodes.
- Carrying case (1).
- Ambulance bracket (EN 1789:2007 certified) (2).
- 4+6 leads patient cable.
- SpO2 probes from Masimo SET.
- EtCO2 probes from Masimo SET.
- “Push Pad” device.
- Rechargeable and non rechargeable battery (3).
- External dual battery charger 12 V (DC) or 220V (AC).
- External battery charger.
- Rechargeable battery for autotraining mode (4).
- “REANIBEX Data Manager, Data Link, Data Cloud” software application to manage and organize all the data gathered.



## Bexen cardio

Osatu, S.Coop.  
Edificio Zearrekobuelta  
Subida de Areitio, 5  
48260 Ermua (Spain)

Tel: +34 943 170 220  
Fax: +34 943 170 227  
info@bexencardio.com  
www.bexencardio.com

