

**Bexen cardio**  
Ready for life

## Reanibex 500

Monitor  
Defibrillator  
Manual/AED with  
Non invasive  
pacing and SpO2



Biphasic Manual and AED  
defibrillator up to 360J

Universal AED algorithm for  
adult and paediatric patients

CPR guidance device in real  
time (Push Pad)

Large colour display (5,7")

Possibility to configure the  
screen on high contrast mode

"Autotraining" mode

# Reanibex 500.

## Rugged and compact.

Pulse Oximetry with Masimo SET technology

Large colour display (5,7")

Disposable electrodes for adult and paediatric patients

Status indicator always visible

CPR guidance with feedback (Push Pad device)

2 USB ports

Electrodes for adult and paediatric patients

3 leads patient cable

## FEATURES

- Universal AED algorithm for adult and paediatric patients.
- Lightweight unit with monitoring capabilities (IEC 60601-2-27).
- 3 leads patient cable.
- Synchronized cardioversion.
- Non-invasive pacing (optional).
- Pulse Oxymetry SpO2 from Masimo SET (optional).
- "Push Pad" device, reports the quality of the frequency and depth of chest compressions in real time (optional).
- Large colour display 5.7".
- Possibility of setting the screen on high contrast mode.
- Configured according to ERC/ AHA guidelines 2015 and its review of 2017.
- Automatic daily selftest ensures the device continual readiness.
- USB connection (2 ports).
- Rugged and durable design for any environment.
- User friendly and very intuitive in use.

## ACCESSORIES

- Adult and paediatric disposable electrodes for defibrillation.
- Adult and paediatric ECG electrodes.
- Carrying case (1).
- Ambulance bracket (EN 1789:2007 certified) (2).
- 3 leads patient cable.
- SpO2 probes from Masimo SET.
- "Push Pad" device.
- Rechargeable and non rechargeable battery (3).
- External battery charger.
- External dual battery charger 12 V (DC) or 220V (AC).
- "REANIBEX Data Manager" Software application to manage and organize all the data gathered.
- Rechargeable battery for autotraining mode (4).



## Bexen cardio

Osatu, S.Coop.  
Edificio Zearrekobuelta  
Subida de Areitio, 5  
48260 Ermua (Spain)

Tel: +34 943 170 220  
Fax: +34 943 170 227  
info@bexencardio.com  
www.bexencardio.com

