

Bexen cardio

Ready for life

Reanibex 700

Monitor Defibrillator Manual



Universal AED algorithm for
adult and paediatric patients

Monitor Defibrillator biphasic
up to 200 J

AED

Non invasive pacing

SpO2

NIBP

Integrated paediatrics paddles

Different
power supply options:

- AC power supply
- Rechargeable battery
- DC power supply

Reanibex 700. Monitor Defibrillator.

Control of the defibrillator and the printer from the paddles

Integrated paediatric paddles

Pulse Oximetry SpO2 from Masimo SET


Non-Invasive Blood Pressure (NIBP) from Suntech Medical

"Compact Flash" memory card

Selector cover: Manual or AED

Non invasive pacing

Visual and acoustic messages



FEATURES

- Universal AED algorithm for adult and paediatric patients.
- Biphasic technology.
- Control of the defibrillator and the printer from the paddles.
- Integrated paediatric paddles.
- Operates from the mains (AC), from a vehicle battery (DC) (optional) and with its internal battery.
- Autotest when switched on and during operation.
- Simple and intuitive.
- AED (optional).
- Non invasive pacing (optional).
- Pulse Oximetry SpO2 from Masimo SET (optional).
- Non-Invasive Blood Pressure from Suntech Medical (optional).

ACCESSORIES

- Disposable electrodes for defibrillation.
- ECG electrodes.
- Carrying case.
- Ambulance bracket (EN 1789:2007 certified) (1).
- 3, 5 or 10 leads patient cable.
- SpO2 probes from Masimo SET.
- Non-Invasive Blood Pressure from Suntech Medical
- Internal paddles (2).
- Rechargeable battery (3).
- "REANIBEX Data Manager" Software application to manage and organize all the data gathered (4).
- External battery charger.
- Connection cable at 12V (DC).



Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobueta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com

