

Bexen cardio

Ready for life

Reanibex 700

Monitor Defibrillator Manual

Monitor Defibrillator biphasic
up to 200 J

Universal AED algorithm for
adult and paediatric patients

Non invasive pacing

SpO₂

NIBP

External adult and paediatric
paddles

50 mm printer


12 leads acquisition

Different
power supply options:

- AC power supply
- Rechargeable battery
- DC power supply



Reanibex 700. Monitor Defibrillator.



50 mm printer

“Compact Flash” memory card

Selector cover:
Manual or AED

Non invasive pacing

Visual and
acoustic messages

Control of the
defibrillator and the
printer from the
paddles

Pulse
Oximetry
SpO2 from
Masimo SET

Integrated paediatric paddles

Non-Invasive Blood Pressure
(NIBP) from Suntech Medical

FEATURES


- AED.
- Universal AED algorithm for adult and paediatric patients.
- Biphasic technology.
- 50 mm printer.
- Operates from the mains (AC), from a vehicle battery (DC) (optional) and with its internal rechargeable battery.
- Control of the defibrillator and the printer from the paddles.
- Integrated paediatrics paddles.
- Autotest when switched on and during operation.
- Simple and intuitive.
- Non invasive pacing (optional).
- Pulse Oximetry SpO2 from Masimo SET (optional).
- Non-Invasive Blood Pressure from Suntech Medical (optional).

ACCESSORIES

- Adult and paediatric disposable pads for defibrillation.
- Adult and paediatric disposable ECG electrodes.
- Carrying case (1).
- Ambulance bracket (EN 1789:2007 + A2:2014 certified) (2).
- 3, 5 or 10 leads patient cable.
- SpO2 probes from Masimo SET.
- Non-Invasive Blood Pressure cuffs from SunTech Medical.
- Rechargeable battery (3).
- External battery charger.
- Connection cable at 12V (DC).

Applications:

- REANIBEX Data Manager PC Application, to store and review ECG data (4).



Bexen cardio

Osatu, S.Coop.
Edificio Zearrekobuelta
Subida de Areitio, 5
48260 Ermua (Spain)

Tel: +34 943 170 220
Fax: +34 943 170 227
info@bexencardio.com
www.bexencardio.com

